

BK5000 LCD Video Biological Microscope



High quality
optical system



Rich optional
accessories



SO, CE, and RoHS
certificate



Novel appearance

Product Details

1. Digital trinocular compound microscope includes a digital imaging system.
2. Professional Trinocular Biological Microscope
3. 9.7" 2048*1536 HD Color Retinal LCD Tablet Camera
4. Tablet camera adopts ANDRROID 5.1.0 operation system
5. The LCD digital biological microscope can connect with PC.
6. Easy and convenient operating procedure save user's time.

Product Parameter

Item Name	BK5000 LCD Video Biological Microscope
Microscope Model	BK5000-TR
Tablet Camera Model	5699
Standard Microscope Outfits	WF10X/20mm eyepiece(dia.23.2mm)
	Infinity Plan 4X/10X/20X/40X/100X(oil) objective
	Inclined 30°, Rotatable 360° trinocular head, Interpupillary Distance:48mm-75mm, Light Distribution: 100: 0 (optional 80:20)
	5 holes nosepiece
	175mm*145mm mechanical stage, travel 78*55mm
	Coaxial coarse and fine focusing knobs on both sides
	Abbe N.A.1.25 condenser with iris diaphragm
	12V/20W halogen lamp,100V-240V
	Blue filter
	0.5X C-mount
9.7 inches LCD tablet camera	Android 5.1.0 operating system
	RK3288 Quad Core 1.8 Hz CPU
	RAM 2GB DDR3
	ROM 8GB
	WI-FI & Bluetooth
	Built- in microscopy software
	Built-in 5mp(2592H*1944px) camera

Product Display



Packing Size

520*480*440mm

Net weight

11 kgs

Gross Weight

13kgs

LCD Tablet Camera

